


















T-Cup®

Urine Drug Test



- CLIA Waived, FDA Approved & Forensic Use Only
- Female User Friendly & Customizable
- Test up to 15 Drugs per Cup

CLIA Waived Drugs

	Amphetamine (AMP)	1000ng/mL 500ng/mL		Methylenedioxymethamphetamine (MDMA)	500ng/mL
	Barbiturates (BAR)	300ng/mL		Morphine (MOP/OPI 300)	300ng/mL
	Benzodiazepines (BZO)	300ng/mL		Opiate (OPI)	2000ng/mL
	Buprenorphine (BUP)	10ng/mL		Oxycodone (OXY)	100ng/mL
	Cocaine (COC)	300ng/mL 150ng/mL		Phencyclidine (PCP)	25ng/mL
	Cannabinoids (THC)	50ng/mL		Propoxyphene (PPX)	300ng/mL
	Methadone (MTD)	300ng/mL		Tricyclic Antidepressants (TCA)	1000ng/mL
	Methamphetamine (MET/mAMP)	1000ng/mL 500ng/mL			

Built-in Adulteration Parameters: Creatinine (CR), Specific Gravity (SG) & pH. Oxidant (OX) & Nitrite (NI) available upon request.

